

Printing date 18.05.2020 Revision: 23.05.2019

### 1 Identification of the substance/mixture and of the company/undertaking

· Product identifier

· Trade name: Opal™ Bond™ MV

· Article number: 71025

· Index number: SDS 368-001.05

· Relevant identified uses of the substance or mixture and uses advised against

Professional Orthodontic Adhesive

· Application of the substance / the mixture Professional Orthodontic Adhesive

· Details of the supplier of the safety data sheet

· Manufacturer/Supplier:

Ultradent Products Inc.

505 W. Ultradent Drive (10200 S)

South Jordan, UT 84095-3942

USA

onlineordersupport@ultradent.com

EC Responsible Person Ultradent Products GmbH

Am Westhover Berg 30

51149 Cologne Germany

Email: infoDE@ultradent.com

Emergency Phone: +49(0)2203-35-92-0

· Further information obtainable from: Customer Service

· Emergency telephone number:

CHEMTREC (NORTH AMERICA) : (800) 424-9300 (INTERNATIONAL) : +(703) 527-3887

### 2 Hazards identification

- · Classification of the substance or mixture
- · Classification according to Regulation (EC) No 1272/2008



GHS07

Skin Sens. 1 H317 May cause an allergic skin reaction.

- · Label elements
- · Labelling according to Regulation (EC) No 1272/2008

The Regulation EC 1272/2008 on classification, labeling and packaging of substances and mixtures (CLP) shall not apply to a medical device in the finished state used in direct physical contact with the human body according to Art. 1.5 (d). Therefore, the product is exempted from the CLP labeling requirements, and no SDS is required by Regulation 1907/2006, Art. 2 (6c), REACH. Therefore, all given data, classification, and information on this SDS are provided solely on a voluntary basis.

- · Hazard pictograms GHS07
- · Signal word Warning
- · Hazard-determining components of labelling:

Diurethane Dimethacrylate

Triethylene Glycol Dimethacrylate

2-Hydroxyethyl Methacrylate

· Hazard statements

H317 May cause an allergic skin reaction.

(Contd. on page 2)

Printing date 18.05.2020 Revision: 23.05.2019

Trade name: Opal™ Bond™ MV

(Contd. of page 1)

#### · Precautionary statements

P101 If medical advice is needed, have product container or label at hand.

P102 Keep out of reach of children.

P103 Read label before use.

*P261* Avoid breathing dust/fume/gas/mist/vapours/spray.

*P272* Contaminated work clothing should not be allowed out of the workplace.

P280 Wear protective gloves.

P302+P352 IF ON SKIN: Wash with plenty of water.

P362+P364 Take off contaminated clothing and wash it before reuse.

P333+P313 If skin irritation or rash occurs: Get medical advice/attention.

*P321* Specific treatment (see on this label).

P501 Dispose of contents/container in accordance with local/regional/national/international regulations.

- · Other hazards
- · Results of PBT and vPvB assessment
- · **PBT:** Not applicable.
- · vPvB: Not applicable.

#### 3 Composition/information on ingredients

- · Chemical characterisation: Mixtures
- · **Description:** Mixture of substances listed below with nonhazardous additions.

· Dangerous components:		
CAS: 72869-86-4	Diurethane Dimethacrylate	>2.5-≤10%
EINECS: 276-957-5	♦ Skin Sens. 1, H317	
CAS: 1344-28-1	Aluminium Oxide	>2.5-≤10%
EINECS: 215-691-6	substance with a Community workplace exposure limit	
CAS: 14808-60-7	Silica Glass	>2.5-≤10%
EINECS: 238-878-4	♠ Acute Tox. 4, H332; Skin Irrit. 2, H315	
CAS: 109-16-0	Triethylene Glycol Dimethacrylate	<i>≤</i> 2.5%
EINECS: 203-652-6	♦ Skin Sens. 1, H317	
CAS: 868-77-9	2-Hydroxyethyl Methacrylate	0-≤2.5%
EINECS: 212-782-2	♦ Skin Irrit. 2, H315; Eye Irrit. 2, H319; Skin Sens. 1, H317	

Additional information: For the wording of the listed hazard phrases refer to section 16.

#### 4 First aid measures

- · Description of first aid measures
- · General information: Immediately remove any clothing soiled by the product.
- · After inhalation:

*This product is a thick paste, therefore inhalation is extremely unlikely.* 

Supply fresh air and to be sure call for a doctor.

*In case of unconsciousness place patient stably in side position for transportation.* 

· After skin contact:

*Immediately wash with water and soap and rinse thoroughly.* 

*If skin irritation continues, consult a doctor.* 

· After eye contact:

Rinse opened eye for several minutes under running water. If symptoms persist, consult a doctor.

- · After swallowing: Do not induce vomiting. Seek medical help if necessary.
- · Information for doctor:
- · Most important symptoms and effects, both acute and delayed No further relevant information available.

(Contd. on page 3)

Printing date 18.05.2020 Revision: 23.05.2019

Trade name: Opal™ Bond™ MV

(Contd. of page 2)

· Indication of any immediate medical attention and special treatment needed No further relevant information available.

#### 5 Firefighting measures

- · Extinguishing media
- · Suitable extinguishing agents:

Foam, dry chemical, carbon dioxide

*Use fire extinguishing methods suitable to surrounding conditions.* 

- · Special hazards arising from the substance or mixture No further relevant information available.
- · Advice for firefighters:
- Protective equipment:

Wear self-contained respiratory protective device.

Wear fully protective suit.

### 6 Accidental release measures

· Personal precautions, protective equipment and emergency procedures

Wear protective equipment. Keep unprotected persons away.

- Environmental precautions: Do not allow to enter sewers/ surface or ground water.
- · Methods and material for containment and cleaning up:

Dispose contaminated material as waste according to item 13.

Ensure adequate ventilation.

· Reference to other sections

See Section 7 for information on safe handling.

See Section 8 for information on personal protection equipment.

See Section 13 for disposal information.

#### 7 Handling and storage

- · Handling:
- Precautions for safe handling: Ensure good ventilation/exhaustion at the workplace.
- · Information about fire and explosion protection: No special measures required.
- · Conditions for safe storage, including any incompatibilities
- Storage:
- Requirements to be met by storerooms and receptacles: Store in a cool location.
- Information about storage in one common storage facility: Not required.
- · Further information about storage conditions:

Protect from exposure to the light.

Protect from heat

See product labelling.

· Specific end use(s) Professional Orthodontic Adhesive

### 8 Exposure controls/personal protection

· Additional information about design of technical facilities: No further data; see item 7.

(Contd. on page 4)

Printing date 18.05.2020 Revision: 23.05.2019

Trade name: Opal™ Bond™ MV

(Contd. of page 3)

#### · Control parameters

#### · Ingredients with limit values that require monitoring at the workplace:

#### 1344-28-1 Aluminium Oxide

WEL (Great Britain) Long-term value: 10\* 4\*\* mg/m³ \*inhalable dust \*\*respirable dust

- Additional information: The lists valid during the making were used as basis.
- · Exposure controls
- · Personal protective equipment:
- · General protective and hygienic measures:

Avoid contact with the eyes and skin.

Immediately remove all soiled and contaminated clothing.

Wash hands before breaks and at the end of work.

Respiratory protection:

In case of brief exposure or low pollution use respiratory filter device. In case of intensive or longer exposure use self-contained respiratory protective device.

· Protection of hands:



#### Protective gloves

The glove material has to be impermeable and resistant to the product/ the substance/ the preparation.

Due to missing tests no recommendation to the glove material can be given for the product/ the preparation/ the chemical mixture.

Selection of the glove material on consideration of the penetration times, rates of diffusion and the degradation

· Material of gloves

The selection of suitable gloves does not only depend on the material, but also on further marks of quality and varies from manufacturer to manufacturer. As the product is a preparation of several substances, the resistance of the glove material can not be calculated in advance and has therefore to be checked prior to the application.

· Penetration time of glove material

The exact breakthrough time has to be found out by the manufacturer of the protective gloves and has to be observed.

- · Eye protection: Not required.
- · **Body protection:** Protective work clothing

#### 9 Physical and chemical properties

- · Information on basic physical and chemical properties
- · General Information
- · Appearance:

Form: Paste
Colour: Whitish
Odour: Acrylic
Odour threshold: Not determined.

· pH-value: Not applicable (non-aqueous)

· Change in condition

Melting point/freezing point: Undetermined.

Initial boiling point and boiling range: Undetermined.

· Flash point: Not applicable.

(Contd. on page 5)

Printing date 18.05.2020 Revision: 23.05.2019

Trade name: Opal™ Bond™ MV

	(Contd. of page
· Flammability (solid, gas):	Not determined.
· Decomposition temperature:	Not determined.
· Auto-ignition temperature:	Product is not selfigniting.
· Explosive properties:	Product does not present an explosion hazard.
· Explosion limits:	
Lower:	Not determined.
Upper:	Not determined.
· Vapour pressure:	Not applicable.
· Density at 20 °C:	2.03 g/cm³
· Relative density	Not determined.
· Vapour density	Not applicable.
· Evaporation rate	Not applicable.
· Solubility in / Miscibility with	
water:	Insoluble.
Partition coefficient: n-octanol/water:	Not determined.
· Viscosity:	
Dynamic:	Not applicable.
Kinematic:	Not applicable.
· Solvent content:	
VOC (EC)	<1 %
Solids content:	<100.0 %
· Other information	No further relevant information available.

## 10 Stability and reactivity

- · Reactivity Stable
- · Chemical stability
- · Thermal decomposition / conditions to be avoided: No decomposition if used according to specifications.
- · Possibility of hazardous reactions: No dangerous reactions known.
- · Conditions to avoid:

Light

Ignition sources

Flames

Heat

- · Incompatible materials: No further relevant information available.
- · Hazardous decomposition products: No dangerous decomposition products known.

### 11 Toxicological information

- · Information on toxicological effects
- · Acute toxicity Based on available data, the classification criteria are not met.

· LD/LC50 values relevant for classification:		
ATE (Acute	Toxicity Estimates)	
Inhalative L	C50/4 h	22.8-37.9 mg/l

(Contd. on page 6)

Printing date 18.05.2020 Revision: 23.05.2019

Trade name: Opal™ Bond™ MV

		(Contd. of page 5
72869-86	-4 Diurethane Dimethaci	rylate
Oral	LD50	>5,000 mg/kg (rat)
1344-28-	l Aluminium Oxide	
Oral	LD50	>5,000 mg/kg (rat)
109-16-0	Triethylene Glycol Dime	thacrylate
Oral	LD50	>5,000 mg/kg (rat)
	LC50 Fish	16.4 mg/l (Fish) (Toxicity to fish)
Dermal	LD50	>2,000 mg/kg (mouse)
868-77-9	2-Hydroxyethyl Methacr	ylate
Oral	LD50	3,275 mg/kg (mouse)
		>5,000 mg/kg (rat)
	LC50 Fish	>100 mg/l (Fish)
Dermal	LD50	>3,000 mg/kg (rabbit)
	LC50(Daphnia magna)	24.1 mg/l (daphnia)

- Primary irritant effect:
- · Skin corrosion/irritation Based on available data, the classification criteria are not met.
- · Serious eye damage/irritation Based on available data, the classification criteria are not met.
- · Respiratory or skin sensitisation

May cause an allergic skin reaction.

- · CMR effects (carcinogenity, mutagenicity and toxicity for reproduction)
- · Germ cell mutagenicity Based on available data, the classification criteria are not met.
- · Carcinogenicity Based on available data, the classification criteria are not met.
- · Reproductive toxicity Based on available data, the classification criteria are not met.
- · STOT-single exposure Based on available data, the classification criteria are not met.
- · STOT-repeated exposure Based on available data, the classification criteria are not met.
- · Aspiration hazard Based on available data, the classification criteria are not met.

#### 12 Ecological information

· Toxicity

· Aquatic toxicity:			
72869-86-4 Diur	72869-86-4 Diurethane Dimethacrylate		
Biodegradability	28 days (Aerobic) (Biodegradability testing)		
EC50	>0.6 mg/l (Algae) (Toxicity to algae)		
	>1.2 mg/l (daphnia) (Toxicity to aquatic invertebrates)		
109-16-0 Triethy	109-16-0 Triethylene Glycol Dimethacrylate		
Biodegradability	28 days (Aerobic) (Biodegradability testing)		
Aqua toxicity	32 mg/l (daphnia) (No Observed Effect Concentration)		
EC50	>100 mg/l (Algae) (Toxicity to algae)		
868-77-9 2-Hydroxyethyl Methacrylate			
EC50	345 mg/l (Algae)		

- · Persistence and degradability No further relevant information available.
- · Behaviour in environmental systems:
- · Bioaccumulative potential No further relevant information available.
- · Mobility in soil No further relevant information available.

(Contd. on page 7)

Printing date 18.05.2020 Revision: 23.05.2019

Trade name: Opal™ Bond™ MV

(Contd. of page 6)

- · Additional ecological information:
- · General notes:

Water hazard class 1 (German Regulation) (Self-assessment): slightly hazardous for water

Do not allow undiluted product or large quantities of it to reach ground water, water course or sewage system.

- · Results of PBT and vPvB assessment
- · **PBT:** Not applicable.
- · vPvB: Not applicable.
- · Other adverse effects No further relevant information available.

### 13 Disposal considerations

- · Waste treatment methods
- · Recommendation Do not allow product to reach sewage system.
- · Uncleaned packaging:

14 Transport information

· Recommendation: Disposal must be made according to official regulations.

· UN-Number · ADR, ADN, IMDG, IATA	not regulated	
· UN proper shipping name · ADR, ADN, IMDG, IATA	not regulated	
· Transport hazard class(es)		
· ADR, ADN, IMDG, IATA · Class	not regulated	
· Packing group · ADR, IMDG, IATA	not regulated	
· Environmental hazards:	Not applicable.	

*Not applicable.* 

Not applicable.

· Transport in	bulk according to Annex II of Marpol	
and the IBC	Code	

· UN "Model Regulation": not regulated

### 15 Regulatory information

· Special precautions for user

· Chemical safety assessment:

Device is biocompatible when used as directed by dental professionals per ISO 10993-1

### 16 Other information

This information is based on our present knowledge. However, this shall not constitute a guarantee for any specific product features and shall not establish a legally valid contractual relationship.

· Relevant phrases

H315 Causes skin irritation.

H317 May cause an allergic skin reaction.

(Contd. on page 8)

Printing date 18.05.2020 Revision: 23.05.2019

Trade name: Opal™ Bond™ MV

(Contd. of page 7)

H319 Causes serious eye irritation. H332 Harmful if inhaled.

- · Department issuing SDS: Regulatory Affairs
- · Contact: Customer Service
- · Abbreviations and acronyms:

ADR: Accord européen sur le transport des marchandises dangereuses par Route (European Agreement concerning the International Carriage of Dangerous Goods by Road)

IMDG: International Maritime Code for Dangerous Goods

IATA: International Air Transport Association

GHS: Globally Harmonised System of Classification and Labelling of Chemicals EINECS: European Inventory of Existing Commercial Chemical Substances

ELINCS: European List of Notified Chemical Substances

CAS: Chemical Abstracts Service (division of the American Chemical Society)

VOC: Volatile Organic Compounds (USA, EU)

LC50: Lethal concentration, 50 percent

LD50: Lethal dose, 50 percent

PBT: Persistent, Bioaccumulative and Toxic

vPvB: very Persistent and very Bioaccumulative

Acute Tox. 4: Acute toxicity - inhalation - Category 4

Skin Irrit. 2: Skin corrosion/irritation – Category 2

Eye Irrit. 2: Serious eye damage/eye irritation – Category 2

Skin Sens. 1: Skin sensitisation – Category 1

GB